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	Interpreting Physicians	Medical Physicists	Radiologic Technologists
Initial Credentialing	<ul style="list-style-type: none"> Medical license, and ABR certification, or AOBR certification, or RCPS certification, or Appropriate training (see text) 	<ul style="list-style-type: none"> ABR certification, or ABMP certification, or State licensure, or State approval 	<ul style="list-style-type: none"> ARRT registered, or ARCRT registered, or State licensure
Initial Training	<ul style="list-style-type: none"> 60 hours Cat I CME (40 hours if qualified before 4/28/99) Eight hours in each modality 	<ul style="list-style-type: none"> Master's in physical science with 20 hours of physics and 20 hours of conducting surveys, or If qualified before 4/28/99, BS in physical science with 10 hours of physics and 40 hours of conducting surveys Eight hours in each modality 	<ul style="list-style-type: none"> If qualified after 4/28/99, 40 hours of training Eight hours in each modality
Initial Experience	<ul style="list-style-type: none"> 240 exams within six months of qualifying date, or If board certified at first opportunity, 240 exams in any six months within last two years of residency 	<ul style="list-style-type: none"> One facility and 10 units, or If qualified before 4/28/99 with BS in physical science, one facility and 20 units 	<ul style="list-style-type: none"> If qualified after 4/28/99, 25 exams under direct supervision
Continuing Experience	960 exams/24 months	Two facilities and six units/24 months	200 exams/24 months
Continuing Education	15 Cat I CME/36 months (six in each modality)	15 CEUs/36 months (include training in each modality)	15 CEUs /36 months (six in each modality)

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Time	General Sessions	Contents
8:30	Introduction and Pre-test	*Viewing Station *Laser Printers *Archiving
9:00	Theory of Digital Mammography	*Current Status of Digital Mammography *Detector Technology *Image Reconstruction and Processing Algorithms *Image Display -Hardcopy and Softcopy *Data Storage -Archiving Issues *Automatic Exposure control *Comparison of Digital to Screen-film Mammography
10:00	Hands-on Acquisition System	*Overview of the Acquisition System *Image Acquisition (AOP vs Manual Techniques etc) *Performing Needle Localizations on the System
11:00	Hands-on Review Workstation	*Buttonology *Soft Copy Review Protocol *Case Examples
12:00	Lunch	
12:30	Quality Control Procedures Overview	
1:00	Hands-on Quality Control Tests	*Technologists Quality Control Tests for Digital Mammography Daily tests *Monitor Cleanliness *Laser Imager Quality Control *Weekly tests *Flat Field Uniformity Test *Phantom Image Quality and Control-to-Noise Ratio Test *Monthly Tests *Visual Checklist *MTF Measurement *AOP(AEC) Mode Check *Quality Tests *Compression Force Test *Repeat Analysis
2:00	Hands-on Group1:Radiologist	*Clinical Case Review *Acceptance Testing of Digital system FFDM Unit Evaluation *Flat Field Uniformity Test Artifact Evaluation *Detector SNR Measurement AOP Mode and SNR check *Phantom Image Quality Test CNR Test *Detector SNR vs Thickness MTF Measurement *Evaluation Assessment System Limiting Resolution *Evaluation of Focal Spot Size Collimation Assessment *kVp Accuracy and Reproducibility Half Value Layer Measurement *Breast Entrance Exposure,Mean Display Monitor Image Quality *Glandular Dose and Radiation
4:30	Summary	
5:00	Closing	Questions , Answers, Post-Test ,and Course Evaluation

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